Requesting Agency/Organization Pinner Community Hospital of Scott Item Details: (Please include all details regarding item specifications and reason for purchase) Sysmex XP300 Hematology Analyzer for Backup and Surse \$10,900 and Date Request to RHC: Date Presented to HCC: HCC Voting Results: Contact Person for this purchase: Title: Laboratory Supervisor Email: Judy Carter Opch Scott.com Phone: 423 286 5350 Shipping Address (if applicable):

One i da TN 37841 *You will need to reference the KET Regional HVA and Gap Analysis to fill in this form*

Item(s) to be Purchased:	Yes	No	Reference/Notes
Please check appropriate category below for this purchase:			
Surge Capacity	V		
Training			
Other	V		Back up - easy to more to alternat
Is this purchase supported by the Regional HVA?			Date of HVA:
Is this purchase supported by lessons learned or best practices from events/exercises?			Explain:
Is there supporting documentation? (AAR, Improvement Plan?)			Attached

18797 Alberta Street

Does the purchase fall within an identified area noted on the KET Coalition Gap Analysis? Y or N

 Healthcare System Pre 	paredness
---	-----------

- 2. Healthcare System Recovery
- 3. Emergency Operations Coordination
- 5. Fatality Management
- 6. Information Sharing
- 10. Medical Surge
- 14. Responder Safety and Health
- 15. Volunteer Management

Function(s) TO BE FILLED OUT BY RHC

- Please send a copy of this to your RHC and you will be given a date to present to KET Coalition meeting.
 To be on the agenda, please submit at least 1 week before the next scheduled meeting.
- Attach sheet with details of purchase (vendor, vendor contact information for payment, cost(quote), etc.)
- Keep in mind, purchases for routine organization operations or to meet joint commission requirements are not an allowable expense.



Quote

DATE 8/26/2015

Quotation # 00698

Customer ID

97661

336 Wilhagan Road

Nashville, TN 37217

Phone 800-626-5275 Fax 615-367-0652

Quotation valid until: 11/24/2015

Prepared by: Robbie Wiedenhoeft

Bill To:

PIONEER HEALTH SERVICES

18797 ALBERTA ST ONEIDA, TN 38741

Ship To: 97661

PIONEER HEALTH SERVICES

18797 ALBERTA ST

ONEIDA, TN 38741

Comments or special instructions

000000000000000000000000000000000000000						romano attinue tituno	et deneración accerciónes
#	Item	Description	UM	Package	Price	Qty	ExtPrice
1	SCAXP300-LAB- KIT	ANALYZER XP300	EA	1 EA / EA	\$10,900.00	1	\$10,900.00


Total: \$10,900.00

Estimated freight is \$150.00.



XP-300™ Automated Hematology Analyzer

Robust and Reliable!



XP-300 - Robust and Reliable

Accurate and reliable

- · Sysmex robustness for the best possible up-time
- · Proven technology for accuracy of results
- · Linearity to support clinical diagnostics and monitoring
- · Sensitive flagging to support diagnosis by the physician

Innovative technology

- Performs rapid and accurate analysis of a 17-parameter CBC
 - Including a 3-part WBC differential with an Absolute Neutrophil Count (ANC)
 - Histograms for RBC, WBC and PLT
 - Automatic floating discriminators
- Produces accurate results comparable to other Sysmex hematology analyzers

Compact and fully integrated

- · Small footprint
- · Fits easily on a laboratory bench
- Ideal as a back-up for Sysmex X-Series[™] and XN-Series[™] analyzers

Easy operation and maintenance

- · Minimal training required
- · Barcoded QC and reagents
- · Auto start-up and shut-down

Safe and secure


- · Non-toxic, biodegradable reagent system
- Two reagents (diluent and lyse) for complete results

Network capability

- Data transfer to LIS via LAN (TCP/IP) or Serial Port (RS-232C)
- · Bi-directional interface



Intuitive touch screen menus simplify operation



QC data is available on board with comprehensive reports and also available throught Insight™ quality control program.

Sysmex Corporation 1-5-1-Wakinohama-Kaigandori Chuo-ku, Kobe 651-0073, Japan Tel. +81 (78) 265-0521 www.sysmex.co.jp Sysmex America, Inc. 577 Aptakisic Road Lincolnshire, IL 60069, U.S.A. Tel. +1 (800) 379-7639 www.sysmex.com/us

Sysmex Canada, Inc. 5045 Orbitor Drive Building 9, Suite 401 Mississauga, ON, L4W 4Y4, Canada Tel. +1 (905) 366-7900 www.sysmex.ca Sysmex Latin America and the Caribbean Rua Joaquim Nabuco 615 - Bairro Cidade Jardim São José dos Pinhais Paraná – Brasil – CEP 83040-210 Tel. +55 (41) 2104-1314 www.sysmex.com.br







©2013 Sysmex America, Inc. MKT-10-1195 07/2013